

HORIZONTAL LAMINAR $\left(\mathbf{N} \right)$ **AIRFLOW CABINETS**

- O Maximum Product/Material Protection
- O Validation Guaranteed Air Cleanliness Class
- O Ergonomics Design / Spacious Work Area




- O Low Noise Level
- O Class 100 / ISO 5 Cleanliness Class



"Laminar Airflow (LAF) cabinet ensures particle-fre eenvironment for product & material/sample protection against airborne dust or microorganisms"

AIRFLOW DIAGRAM



Metisafe® Horizontal LAF cabinets works under air filtration and indoor exhausting principle. Room air pass through pre-filters that blocks coarse particles at devize top and reaches to work cabinet area with air velocity of 0,40 m/s $(\pm\% 5)$ by passing through HEPA filter. Built in laminator in fan filter module supplies the airflow at homogenous and laminated way. Designs of working cabinet and air diffusion assures trimmed airflow and sweeping dead zones. Air moves in one way in working cabinet by exhausted to room at front by preventing room air move into the working zone. The automatic air compensation through microprocessor control system keeps steady airflow rate at set point.

Air Passing through the Prefilter



- O Homogenous and cruise control airflow
- O Advanced filter compensation system keeping steady air velocity under increased filter resistance
- O Negative plenum Fan Filter Unit suitable for non-stop operation
- O Homogenous illumination
- O Timer controlled effective disinfection by UV Lamp
- O D.O.P. test ports

Main Body and Design Features

- O Main body from antibacterial electrostatic paint coated on steel
- O <7° angle back sloped front panel for work comfort and increased working chamber vision
- O Easy reachable electronic control & electrical part located inside at the front panel isolated from the work chamber.
- O UV resistant glass




MICROPROCESSOR CONTROL UNIT


Color Screen Touch Control Panel: Password protected user interface, Device on/off button, Standby/Normal modes operation indicator, Motorized front window level adjustment and one-button automatic front window level adjustment, UV and illumination lamp On/Off buttons, Automatic switch to Standby mode and UV lamp operation when the windshield is closed, User interface password change, Password protected technical service maintenance/calibration, Date and time information

Audio-Visual Warning System: Audible visual alarm button and temporary mute, Improper front window position warning, Improper air velocity warning, Filter, pre-filter and UV lamp replacement warning, Technical service maintenance warning


Large Information Display: Downflow air velocity rate, Fan working capacity, Fan working diagram, UV lamp timer setting, Cabinet, UV and illumination running times




Easy access, isolated control unit from the contaminated are



Monoblock with rounded corners stainless steel work plate



Control panel main display screen



Plenum cast design creates negative pressure around contaminated work area and ensures zero leakage



OPTIONALS

MOP.ADI	Lux Adjustable Illumination
MOP. FWS	Front Window Motion Safety System
MOP.HTS	Humidity/Temperature Sensor
MOP.RES	Remote Access
MOP.CFL.EU4	Active Carbon Impregnated Pre-filter
MOP.FLT.ULP	ULPA Filter
MOP.MIE	Microscope Entegration
MOP.PMS	Particle Measure System
MOP.MDS	Manual Decontamination System
MOP.WCS	Working Cabinet Suspension Apparatus
MOP.BWT	Precision Balance with Vibration
	Damping Weighing Table

[+'

MAC MAC MAC MAC MAC MAC MAC MAC MAC

MAC MAC MAC MAC

·#	ACCESSORIES
.HAS	Height Adjustable Support Stand
.WSS	Wheeled Support Stand
.UVC	Portable Germicidal UV Lamp
.MIS	Microbiological Air Sampling Apparatus
.ASA	Arm Support Apparatus
UPS.	Uninrerruptible Power Supply (UPS)
.ADH	Acrobat Document Holder
AMH.	Acrobat Monitor Holder
.INC	Incinerator
DID	Disinfection Dispenser
ECM.	Document/Equipment Cleaning Module
.VLF	Service Valves
.FOS	Foot Stand



Apparatus (IV)

· · ·

Incinerator



Arm Support Apparat








Height Adjustable


Support Stand



Microscope Integration

METİSAFE HORIZONTAL TYPE LAF CABINET TECHNICAL DRAWING

0 0



Model			LFH-90	LFH-120	LFH-150	LFH-180	
DIMENSI	ONS AND WEIG	іНТ					
Internal Dimensions (E x D x Y) mm External Dimensions (E x D x Y) mm Support Stand (Ex Dx Y) mm Working Table Height (mm) Cabinet Weight		•	875 x 550 x 586	1180 x 550 x 586	1490 x 550 x 586	1790 x 550 x 586	
		995 x 865 x 1345	1300 x 865 x 1345	1610 x 865 x 1345	1910 x 865 x 1345		
		995 x 825 x 720	1300 x 825 x 720	1610 x 825 x 720	1910 x 825 x 720		
		750 mm					
			150 kg	186 kg	227 kg	272 kg	
	Support Stand Weight		30 kg	34 kg	38 kg	42 kg	
-	Packaged Total Weight (With Support Stand)		215 kg	260 kg	310 kg	364 kg	
	Pack/Palette Dimensions (With Support Stand) (WxLxH) mm		1140 x 1105 x 1780	1450 x 1105 x 1780	1760 x 1105 x 1780	2070 x 1105 x 178	
HEPA FILTER and AIRFLOW DYNAMICS							
Airflow \	/elocitv	tity Down Flow 0.45 m/s ± %20					
		Pre-Filter	EU4				
Filters (E	Filters (EN 1822)	Main Filter	H14	H14 HEPA, 0,3 µm paticle %99.995 < filtration efficiency			
Airflow	Rate m³/h		615	770	975	1165	
	-	EN ISO 14644-3	015			1105	
Working Cleaning		US FED 209E	< ISO 5 < Class 100				
0.04	0.000	031202032					
ENERGY							
ENERGY		Fan Filter Unit	90 W	120 W	136 W	148 W	
		Illumination	18 W	36 W	36 W	36 W	
	Energy Consumption 230 VAC- 50Hz	UV Lamp	18 W	30 W	36 W	36 W	
230 VAC	- 50HZ	Electrical Socket					
		Power / Current	2x1150 W / 2x5 A				
	wer / Current		2426 W/10,5 A	2486 W/10,8 A	2508 W/10,9 A	2520 W/10,95 A	
PowerC	Power Consumption Under Normal Operating Conditions (Illumination and Fan Motor)		106 W	152 W	168 W	180 W	
		·					
ERGONO	MICS and COM	FORT					
Front Wi	indow Aperture	(Fixed) - Optional		40	00 mm		
Front Wi	indow	Standard-Preset	200 ± 10mm				
Aperture	e - Motorized	Max. Height	410 ± 10mm				
Noise Le	vel	Normal Mod	≤ 60 dB				
Illumination			≥ 750 Lux				
t t							
WORKIN	IG AREA						
		Standard	304 Stainless Steel				
Working	Table	Optional	316 Stainless Steel				
	Front Window Thickness			≥ 6 mm			
	indow Thickness			~	6 mm		

Working Area

CE, ISO9001-2015 Quality Management Certificate, ISO14644 Cleanroom Class Compatibility, EN61010 Electric and Electronic Safety Compatibility, Accredited Testing Institution Validation Guarantee in case the installation is opareted by Metisafe® Certified Technical Personnel



Telephone : +90(312) 397 64 99 | Fax : +90(312) 395 30 77 Saray Mah. 185.Cad. No:4 Kahramankazan, 06980 ANKARA-TURKEY www.metisbio.com / info@metisbio.com

METISAFE® is trademark of Metis Biotechnology

This content is for informational purposes only and all rights are reserved. Metisafe® has the right for any kind of content and information updates.